How is a pH meter calibrated?

Mr. M. Gonzalez

Do now (5 minutes)

► Grab an Apron.

There are some that are cut short for our classmates. If you can use a longer apron, then use it and leave the ones that are cut off for others.

THANK YOU!

Do Now: Germination Rate Calculation

Today is day 6 from when we began germination.

Calculate today's "real rate."

Day	Sprouts	Sponge Total	Hyp. Rate %	5% error	Real Rate %	Support Y / N
6		32				
7		32				
8		32				
9		32				
10		32				

How dobews is yamphilerbie telli nahiabbattéon?

Objectives

I can calibrate a pH meter.I can understand what instrument error is.

Materials

Calibration

- Hanna Combo Meter (pH meter)
- pH buffer capsules for pH 4, 7, and 10
- Distilled water
- Graduated cylinder
- Stirring rod optional
- Jar with cap (greater than 100mL)
- Colored labeling tape
- Paper towels
- Apron and goggles
- Non latex gloves





Fig. 03. 1000mL beaker and squirt bottle.

Fig. 04. Probe storage solution and probe.

Combo Waterproof

Calibrating Meter

- Take the cap off
- Rinse off the probe
- Turn meter on
- Press "Set" button until reading pH
- Press and hold "Mode" button until it reads "Off" turns to "Cal"
- Place probe into standard 7.





Calibrating Meter

- When meter changes to 4, rinse off probe and place in the 4 standard
- Wait until the "Cal" display stops blinking
- Meter is now calibrated and ready for use.

Activity

Create the following data table:

Meter	pH
Yours	
Other Group 01	
Other Group 02	
Other Group 03	

Activity

Measure the pH of the water sample with your meter.

Add your reading to the table.
Get another three readings from other groups.
Explain why the instruments are not all reading the same pH.

Measuring pH

- Turn on the meter by pressing the Mode button firmly.
- If the letters "pH" do not appear at the top righthand portion of the screen press the set button firmly and continuously until it does.





Activity

Measure the pH of the water sample with your meter.

Add your reading to the table.
Get another three readings from other groups.
Explain why the instruments are not all reading the same pH.

Measuring pH

- If necessary, calibrate the meter.*
- Place the probe into the solution. Wait until the primary display stabilizes (15s) before you read the value.



Measuring pH

- If it doesn't stabilize, look at the value, look away, count to three, repeat twice and average out the value.
- Unless the value keeps on descending or ascending and the value looks intuitively way off. In this case make a note of the problem and tell the instructor.



Measuring pH

- Rinse the probe after each sample solution.
- Read other assigned samples.
- After your final rinse place 4 drops of pH storage solution into the little cup located inside the probe cap.





Clean-up

- Recap the probe FIRMLY by lining up the lettering on the cap and the face of the meter.
- Wait for clean-up instructions.
- Clean up your area, dump your rinse water out, dispose of gloves, and return materials.
- Begin work on your bar graph and analysis questions.

Instrument Error

Why do you think different meters gave different pH values for the same solution?
Take notes on the answers given in class about instrument error.

Do now

Read "Measuring the World" by Chang (2005).
Be prepared to answer what the standard for a meter is?

Also, how does this knowledge of how we define what we measure change your understanding of measurement in general.